CONTRASTE DE FASES EN TORRETA TURRET PHASE CONTRAST EQUIPMENT

REFERENCIA 90137022

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used, the 10x position of the turret must match the white spot at the front).

◆ Change one of the eyepieces of the microscope for the centring eyepiece supplied with the phase contrast equipment.

◆ By turning the centring eyepiece, focus sharply the bright ring and the phase contrast plate (dark ring).

◆ Proceed to centre the bright ring over the phase contrast plate. For this, move gently the corresponding diaphragm located under the turret; by observing through the centring eyepiece check that the bright ring coincides with the dark ring (*figure 4*). Focus both rings with the centring eyepiece.

◆ Remove the centring eyepiece and replace the eyepiece that is going to be used. The microscope is ready for observation.

Note: if the image is not sharp, maybe the sample is too diffusing and is not possible to be observed with a phase contrast system.

REMEMBER

For a phase contrast observation, the different structures must be unstained or slightly stained; they must be extremely thin and must not present thick walls.

The optical surfaces of the microscope as well as the slide must be perfectly clean.

The light intensity must be higher than for bright field observations due to the ring-shaped reduction of the incident light and above all, to the absorption of the phase contrast plate.


Figure 3



Figure 4

1. STANDARD CONFIGURATION

Please check that all the accessories are enclosed with the equipment:

- Phase contrast turret with condenser
- Phase contrast plan achromatic objectives: 10x, 20x, 40x (R) y 100x (R) (I)
- Centring eyepiece
- Screw driver
- Aluminium case

2. INSTALLATION

◆ Prepare the nosepiece (1.1) of the microscope with the phase contrast and usual objectives that are going to be used.

◆ Remove the condenser and diaphragm (1.2) after loosening the holding screw.

◆ Insert the condenser that is joined to the phase contrast turret; tighten well the screw with the screw driver supplied since the turret is quite heavy and could fall (figure 2).

◆ Put a slide on the microscope stage and focus the image as usual; for this purpose, turn the phase contrast turret until the spot marked as "B" matches the white spot at the front of the turret. To focus the image use the 4x or 10x objective that is normally used with the microscope.

◆ For a phase contrast observation of the sample, put the desired phase contrast objective at the central position (visualization position) of the nosepiece. Select the position of the turret (figure 3) that matches the magnification of the selected objective (for example, if the 10x objective is going

Figure 1



Figure 2

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
Figura 3


Figura 4

Thank you for choosing this equipment. We sincerely wish that you enjoy your Zuzi turret phase contrast equipment. We highly recommend looking after this equipment according to what is stated in this manual.

Zuzi develops its products according to the CE marking regulations as well as emphasizing the ergonomics and security for its user.

The correct using of the equipment and its good quality will permit you to enjoy this equipment for years.

The improper use of the equipment can cause accidents and electric discharges, circuit breakers, fires, damages, etc. Please read the point of Maintenance, where we expose the security notes.

TO GET THE BEST RESULTS AND A HIGHER DURATION OF THE EQUIPMENT IT IS ADVISABLE TO READ THOROUGHLY THIS MANUAL BEFORE OPERATING WITH THE EQUIPMENT.

Please bear in mind the following:

- ♦ This manual is inseparable from the Zuzi turret phase contrast equipment, so it should be available for all the users of this equipment.
- ♦ You should carefully handle the equipment avoiding sudden movements, knocks, free fall of heavy / sharp objects on it. Avoid spilling liquids inside the equipment.
- ♦ Never dismantle the different pieces of the equipment to repair it yourself, since it could produce a defective use of the whole equipment and a loss of the product warranty, as well as injuries on people that handle the equipment.
- ♦ If you have any doubt about setting up, installation or functioning do not hesitate in contacting your wholesaler. You can also tell us any doubts or suggestions you have by contacting Zuzi Technical Assistance Department by email to asistencia@auxilab.es or by telephone: +34 807 117 040 (0.30 Euros/min).
- ♦ This equipment is protected under the Warranties and consumer goods regulation (10/2003).
- ♦ Overhaul is not covered by the equipment warranty.
- ♦ Operations made by non-qualified staff will automatically produce a loss of the warranty.
- ♦ Please make sure you keep the invoice, either for having the right to claim or asking for warranty coverage. In case you have to send the equipment to Zuzi Technical Assistance Department you should enclose the original invoice or a copy as guarantee.
- ♦ Manufacturer reserves the right to modify or improve the manual or equipment.

ATTENTION!! IF EQUIPMENTS ARE NOT PROPERLY CLEAN AND DISINFECTED THEY WOULD NOT BE ALLOWED TO REPAIR BY OUR TECHNICAL SERVICE.